


# TD4A-2227-UK MKII Mini Telephone Monitor

Part of the SA3000 Paging System for deaf and deaf-blind people



## Getting Started


1. Install a 9 volt PP3 Alkaline or Lithium battery. **ONLY USE ALKALINE OR LITHIUM TYPE BATTERIES.**
2. When a battery is first fitted the indicator light flashes red then green several times, there is a brief pause and it then flashes green several times. The number of green flashes in this second sequence indicates the voltage in the battery. If the light flashes red then the battery is flat and should be replaced. To repeat the battery test, remove the battery then press the Push Button. When the battery is reconnected the battery will be retested automatically. Automatic low battery signals are sent to the Pager or SignWave when the battery runs low. In addition, the indicator light will shine red when the unit is activated or when the push button is pressed, if the battery is low.
3. Plug the Mini Telephone Monitor into the telephone socket. It may be necessary to use a double adaptor if a spare socket is unavailable. If connecting to a broadband line **A MICRO FILTER MUST BE CONNECTED TO THE LINE BEFORE CONNECTING THE TELEPHONE MONITOR**
4. Press the push button, the green light on the front on the monitor will light for a number of seconds; check that a signal is received by the Pager or SignWave. Then ring the telephone line to ensure the monitor activates.
5. It is possible to use the Telephone Monitor to detect incoming mobile calls and text messages. It is possible to do this and detect incoming landline calls using a single Telephone Monitor. Control switch 6 needs to be moved to the On position and the mobile phone placed close to the Telephone Monitor. It may be necessary to experiment to find the best place to locate the mobile in relation to the monitor. A good starting point is adjacent to the monitor as shown below. The mobile call detector in the monitor is located at the bottom of the unit near to the cable exit.
6. It is possible to alter the symbol that lights up on the Pager & SignWave when the Telephone Monitor is activated. If you have more than one telephone line you will require two monitors. If key 5 is moved to the on position on one monitor a person to person signal will be sent instead of the usual telephone signal allowing you to determine which line is ringing.

### Control Switch settings

7. Control Keys 1-4 are for system coding. These must correspond to the settings on the Pager or SignWave.
8. Control Key 5 will cause the Person to Person symbol to light on the Pager & SignWave.
9. Control Key 6 is used to alert of incoming mobile calls & text messages.
10. Control Key 7 not used.
11. Control Key 8 not used.



ON



9 100% Recycled Paper

**New 2 Year Warranty**

On All Silent Alert Products produced after 8<sup>th</sup> June 2009

**NOT SUITABLE FOR CHILDREN UNDER 36 MONTHS BECAUSE OF SMALL PARTS. A MICRO FILTER MUST BE CONNECTED TO THE LINE BEFORE CONNECTING THE TELEPHONE MONITOR.**

The SA3000 System comprises a SignWave Receiver (with Base) or Pager Unit (with either an Alarm Clock Charger or Trickle Charger), a Pillow Pad and at least one transmitter option. Clofield Ltd., Chesterfield, England.

Help Line 0800-387-397, Txt 01246-572141, Fax 01246-450789 e-mail [deaf.info@silent-alert.co.uk](mailto:deaf.info@silent-alert.co.uk)