

CCT3A-2215-IN Trickle Charger

Part of the SA3000 Care Call System



Getting Started

1. Locate the Trickle Charger on a bedside table or desk top and connect the Power Supply
2. Plug the Power Supply into a 13 amp mains socket. Leave the Charger switched on at all times as it includes an emergency back-up battery.
3. Ensure the Pager unit is switched on and dock it onto the charging plug ensuring the Pager is pushed fully down. Charge the Pager for 24 hours.
4. Install the Care Call Monitor units (Part Nos CCM3A-2214-EU, CCM3A-2216-EU-PIR or CCF3A-2219-EU) according to their instructions.
5. Whilst the Pager is in the charger slot operate each Care Call Monitor in turn. The Pager will respond by lighting the relevant key and activating the internal sounder.
6. To reset the Pager, press and hold the key that is lit. Each signal will automatically cancel after 1-2 minutes.
7. The battery key on the Pager will not flash whilst it is in the Charger.
8. The Pager must be recharged each night and can be left on charge permanently.

DO NOT IMMERSE THE TRICKLE CHARGER IN WATER OR ANY OTHER LIQUID. DO NOT SUBJECT THE UNIT TO EXCESSIVE MECHANICAL SHOCK. THE TRICKLE CHARGER SHOULD BE PLACED ON A SMOOTH HARD SURFACE AND AWAY FROM OTHER OBJECTS TO ENSURE CORRECT VENTILATION. DO NOT COVER UNIT OR PLACE ON SOFT OR CARPETED SURFACES.