


## CL4A-2212-EU Mini Chair Leaving Monitor

Part of the Silent Alert Paging System for deaf and deaf-blind people



### Getting Started

1. Install a 9 volt PP3 Alkaline or Lithium battery. **ONLY USE ALKALINE OR LITHIUM TYPE BATTERIES.**
2. When a battery is first fitted the indicator light flashes red then green several times, there is a brief pause and it then flashes green several times. The number of green flashes in this second sequence indicates the voltage in the battery. If the light flashes red then the battery is flat and should be replaced. To repeat the battery test, remove the battery then press the Push Button. When the battery is reconnected the battery will be retested automatically. Automatic low battery signals are sent to the Pager or SignWave when the battery runs low. If the battery is low the indicator light will shine red when the unit is activated by the Pad or after the push button is pressed and released.
3. To test the unit press then release the Push Button. There will be a short delay then the indicator light will shine green and a signal is sent.
4. The unit will by default cause the Person to Person symbol light on the Pager or SignWave. Key 5 can be used to change the signal to light the Burglar symbol.
5. Plug a Care Call Chair Leaving Pad (PM5) to the socket on the Monitor. Place the Pad and Monitor under the cushion of a chair, ensuring that the pad is centrally located and the monitor is placed towards the back of the chair. It is important that the Monitor push button is not depressed when the user is seated on the chair. As the user rises from the chair the Monitor will send a radio signal to the Pager or SignWave.
6. It is possible to increase the range of the system. To do this move key 5 on the Pager or SignWave to the on position. In the case of the Pager it will be necessary to switch the Pager off for five seconds before moving key 5 to the ON position and then switch the Pager back on again.

#### Control Switch settings (move away from the numbers on the plastic case to switch On)

1. Control Keys 1-4 are for system coding. These must correspond to the settings on the Pager or SignWave.
2. Control Key 5 change to illuminate the Burglar Symbol
3. Control Key 6 System channel override
4. Control Key 7 **must be left in the On position**
5. Control Key 8 Not used

### New 2 Year Warranty

On All Silent Alert Products produced after 8<sup>th</sup> June 2009

9 100% Recycled Paper

**NOT SUITABLE FOR CHILDREN UNDER 36 MONTHS BECAUSE OF SMALL PARTS. DO NOT APPLY MAINS POWER OR ANY VOLTAGE IN EXCESS OF 30 VOLTS AS DAMAGE WILL RESULT.**